

LACTOSCAN DISPLAY MODULE

Clear Measurement Display • Physical Button Control • Reliable Operation



**NON-TOUCH
LCD DISPLAY**

The display module is used to show measurement results and analyzer status, while physical buttons are used for selecting/starting functions and navigation on supported Lactoscan configurations.

Lactoscan systems can display measurement parameters such as:

- Fat
- SNF
- Density
- Protein
- Lactose
- Water / Added Water
- Temperature

Additional parameters may be available depending on model and configuration.




Easy-to-read results

Get clear and accurate measurement data at a glance.



Dedicated control buttons

Quick and convenient navigation and function selection.



Rear View (PCB)



Clear Readout



Physical Key Control



Model-Compatible
Operation